


ENDOVASCULAR AORTIC CONFIGURATIONS

Together, improving life



GORE® EXCLUDER® AAA Endoprosthesis with C3® Delivery System

Trunk-Ipsilateral Leg Endoprosthesis

Catalogue number	Intended aortic inner diameter (mm)	Aortic endoprosthesis diameter (mm)	Length to flow divider (cm)	Intended iliac inner diameter (mm)	Iliac endoprosthesis diameter (mm)	Overall device length (cm)	Recommended sheath size (Fr)
RLT231212	19–21	23	4	10–11	12	12	16
RLT231214	19–21	23	4	10–11	12	14	16
RLT231216	19–21	23	4	10–11	12	16	16
RLT231218	19–21	23	4	10–11	12	18	16
RLT231412	19–21	23	4	12–13.5	14.5	12	16
RLT231414	19–21	23	4	12–13.5	14.5	14	16
RLT231416	19–21	23	4	12–13.5	14.5	16	16
RLT231418	19–21	23	4	12–13.5	14.5	18	16
RLT261212	22–23	26	4	10–11	12	12	16
RLT261214	22–23	26	4	10–11	12	14	16
RLT261216	22–23	26	4	10–11	12	16	16
RLT261218	22–23	26	4	10–11	12	18	16
RLT261412	22–23	26	4	12–13.5	14.5	12	16
RLT261414	22–23	26	4	12–13.5	14.5	14	16
RLT261416	22–23	26	4	12–13.5	14.5	16	16
RLT261418	22–23	26	4	12–13.5	14.5	18	16

GORE® EXCLUDER® AAA Endoprosthesis with C3® Delivery System

Trunk-Ipsilateral Leg Endoprosthesis (continued)

Catalogue number	Intended aortic inner diameter (mm)	Aortic endoprosthesis diameter (mm)	Length to flow divider (cm)	Intended iliac inner diameter (mm)	Iliac endoprosthesis diameter (mm)	Overall device length (cm)	Recommended sheath size (Fr)
RLT281212	24–26	28.5	4	10–11	12	12	18
RLT281214	24–26	28.5	4	10–11	12	14	18
RLT281216	24–26	28.5	4	10–11	12	16	18
RLT281218	24–26	28.5	4	10–11	12	18	18
RLT281412	24–26	28.5	4	12–13.5	14.5	12	18
RLT281414	24–26	28.5	4	12–13.5	14.5	14	18
RLT281416	24–26	28.5	4	12–13.5	14.5	16	18
RLT281418	24–26	28.5	4	12–13.5	14.5	18	18
RLT311413	27–29	31	5	12–13.5	14.5	13	18
RLT311415	27–29	31	5	12–13.5	14.5	15	18
RLT311417	27–29	31	5	12–13.5	14.5	17	18
RLT351414	30–32	35	6	12–13.5	14.5	14	18
RLT351416	30–32	35	6	12–13.5	14.5	16	18
RLT351418	30–32	35	6	12–13.5	14.5	18	18


GORE® EXCLUDER® AAA Endoprosthesis with C3® Delivery System

Aortic Extender Endoprosthesis

Catalogue number	Intended aortic vessel diameter (mm)	Endoprosthesis diameter (mm)	Endoprosthesis length (cm)	Recommended introducer sheath (Fr)
PLA230300	19–21	23	3.3	16
PLA260300	22–23	26	3.3	16
PLA280300	24–26	28.5	3.3	16
PLA320400	27–29	32	4.5	17
PLA360400	30–32	36	4.5	18

GORE® EXCLUDER® AAA Endoprosthesis with C3® Delivery System

Aortic Extender Endoprosthesis and Trunk-Ipsilateral Leg Endoprosthesis



GORE® EXCLUDER® Conformable AAA Endoprosthesis with ACTIVE CONTROL System

Trunk-Ipsilateral Leg Endoprosthesis


Catalogue number	Intended aortic vessel diameter (mm)	Aortic endoprosthesis diameter (mm)	Length to flow divider (cm)	Intended iliac vessel diameter (mm)	Iliac endoprosthesis diameter (mm)	Overall device length (cm)	Recommended sheath size (Fr)
CXT201212E	16–18	20	5.5	10–11	12	12	15
CXT201412E	16–18	20	5.5	12–13.5	14.5	12	15
CXT231212E	19–21	23	5.5	10–11	12	12	15
CXT231412E	19–21	23	5.5	12–13.5	14.5	12	15
CXT261212E	22–23	26	5.5	10–11	12	12	16
CXT261412E	22–23	26	5.5	12–13.5	14.5	12	16
CXT281212E	24–26	28.5	5.5	10–11	12	12	16
CXT281412E	24–26	28.5	5.5	12–13.5	14.5	12	16
CXT321414E	27–29	32	6.5	12–13.5	14.5	14	18
CXT361414E	30–32	36	6.5	12–13.5	14.5	14	18

Aortic Extender Endoprosthesis

Catalogue number	Intended aortic vessel diameter (mm)	Endoprosthesis diameter (mm)	Endoprosthesis length (cm)	Recommended introducer sheath (Fr)
CXA200005E	16–18	20	4.5	15
CXA230005E	19–21	23	4.5	15
CXA260005E	22–23	26	4.5	15
CXA280005E	24–26	28.5	4.5	16
CXA320005E	27–29	32	4.5	18
CXA360005E	30–32	36	4.5	18

GORE® EXCLUDER® Conformable AAA Endoprosthesis with ACTIVE CONTROL System

Aortic Extender Endoprosthesis and Trunk-Ipsilateral Leg Endoprosthesis



GORE® EXCLUDER® AAA Endoprosthesis with C3® Delivery System

GORE® EXCLUDER® Conformable AAA Endoprosthesis with ACTIVE CONTROL System

Contralateral Leg Endoprosthesis

Catalogue number	Intended iliac vessel diameter (mm)	Iliac endoprosthesis diameter (mm)	Contralateral leg length (cm)	Recommended introducer sheath (Fr)
PLC121000	10–11	12	10	12
PLC121200	10–11	12	12	12
PLC121400	10–11	12	14	12
PLC141000	12–13.5	14.5	10	12
PLC141200	12–13.5	14.5	12	12
PLC141400	12–13.5	14.5	14	12
PLC161000*	13.5–14.5	16	9.5	12
PLC161200*	13.5–14.5	16	11.5	12
PLC161400*	13.5–14.5	16	13.5	12
PLC181000*	14.5–16.5	18	9.5	12
PLC181200*	14.5–16.5	18	11.5	12
PLC181400*	14.5–16.5	18	13.5	12
PLC201000*	16.5–18.5	20	9.5	12
PLC201200*	16.5–18.5	20	11.5	12
PLC201400*	16.5–18.5	20	13.5	12
PLC231000*	18.5–21.5	23	10	14
PLC231200*	18.5–21.5	23	12	14
PLC231400*	18.5–21.5	23	14	14
PLC271000*	21.5–25	27	10	15
PLC271200*	21.5–25	27	12	15
PLC271400*	21.5–25	27	14	15

* PLEASE NOTE: All 16, 18, 20, 23, 27 mm diameter Contralateral Leg Endoprotheses can be used as Iliac Extenders.


GORE® EXCLUDER® AAA Endoprosthesis with C3® Delivery System

GORE® EXCLUDER® Conformable AAA Endoprosthesis with ACTIVE CONTROL System


Iliac Extender Endoprosthesis

Catalogue number	Intended iliac vessel diameter (mm)	Iliac endoprosthesis diameter (mm)	Contralateral leg length (cm)	Recommended introducer sheath (Fr)
PLL161007	8–9	10	7	12
PLL161207	10–11	12	7	12
PLL161407	12–13.5	14.5	7	12


Contralateral Leg Endoprosthesis



Iliac Extender Endoprosthesis



- All large diameter Contralateral Leg Endoprostheses (16, 18, 20, 23, 27 mm) can be used as Iliac Extenders.
- The 27 mm x 10 cm Contralateral Leg Endoprostheses does not have a radiopaque marker 6 cm above the distal end.
- Traditional sizes of Iliac Extender (10, 12, 14.5 mm) are not compatible with large diameter (16, 18, 20, 23, 27 mm) Contralateral Leg Endoprostheses.



GORE® EXCLUDER® Iliac Branch Endoprosthesis


Iliac Branch Component

Catalogue number	Intended external iliac vessel diameter (mm)	Diameter iliac end (mm)	Iliac branch component length (cm)	Recommended sheath size (Fr)
CEB231010	6.5–9	10	10	16
CEB231210	10–11	12	10	16
CEB231410	12–13.5	14.5	10	16

Internal Iliac Component

Catalogue number	Intended internal iliac vessel diameter (mm)	Diameter iliac end (mm)	Internal iliac component length (cm)	Recommended sheath size (Fr x cm)
HGB161007	6.5–9	10	7	12 x 45
HGB161207	10–11	12	7	12 x 45
HGB161407	12–13.5	14.5	7	12 x 45

GORE® EXCLUDER® Iliac Branch Endoprosthesis




GORE® EXCLUDER® Iliac Branch Endoprosthesis


Deployment sequence




1 Establish sheaths and up-and-over throughwires via snare technique.




2 Advance Iliac Branch Component over both wires on GORE® EXCLUDER® Iliac Branch Endoprosthesis treatment side.



3 Withdraw sheath below GORE® EXCLUDER® Iliac Branch Endoprosthesis and align radiopaque marker ring above internal iliac ostium.




4 Deploy Iliac Branch Component (pull white outer deployment knob) and advance 12 Fr x 45 cm flexible sheath over bifurcation and into internal iliac gate.




5 Cannulate internal iliac artery with third guidewire. Partially deployed Iliac Branch Component can be repositioned by: A) Rotating up to 90° in either direction, or B) Moving distally, to facilitate access to the internal iliac artery.

GORE® EXCLUDER® Iliac Branch Endoprosthesis


Deployment sequence




6 Advance Internal Iliac Component and align short Internal Iliac Component radiopaque marker with long Iliac Branch Component radiopaque marker.




7 Deploy Internal Iliac Component and balloon overlap with Iliac Branch Component and seat within internal iliac artery. Remove up-and-over throughwire.



8 Deploy external iliac leg of Iliac Branch Component (pull gray inner deployment knob).



9 Exchange the 12 Fr flexible sheath for an 18 Fr sheath. Advance and deploy the Trunk-Ipsilateral Leg Endoprosthesis of the GORE® EXCLUDER® Device.



10 Advance and deploy 23 or 27 mm diameter Contralateral Leg as bridging component.



GORE® TAG® Conformable Thoracic Stent Graft with ACTIVE CONTROL System

Catalogue number	Intended aortic diameter (mm)	Proximal diameter (mm)	Distal diameter (mm)	Endoprosthesis length (cm)	Oversizing range (%)	Partially uncovered stent length (mm)	Recommended sheath size (Fr)
TGM212110E	16–19.5	21	21	10	8–31	3	18
TGM262110E	19.5–24 / 16–19.5	26	26	10	8–33	4	20
TGM262610E	19.5–24	26	26	10	8–33	4	20
TGM282810E	22–26	28	28	10	8–27	4	20
TGM282815E	22–26	28	28	15	8–27	4	20
TGMR312610E	24–29 / 19.5–24	31	31	10	7–33	4	20
TGMR313110E	24–29	31	31	10	7–29	4	20
TGMR313115E	24–29	31	31	15	7–29	4	20
TGMR313120E	24–29	31	31	20	7–29	4	20
TGM343410E	27–32	34	34	10	6–26	5	22
TGM343415E	27–32	34	34	15	6–26	5	22
TGM343420E	27–32	34	34	20	6–26	5	22
TGMR373710E	29–34	37	37	10	9–28	5	22
TGMR373715E	29–34	37	37	15	9–28	5	22
TGMR373720E	29–34	37	37	20	9–28	5	22
TGMR404010E	31–37	40	40	10	8–29	6	22
TGMR404015E	31–37	40	40	15	8–29	6	22
TGMR404020E	31–37	40	40	20	8–29	6	22
TGM454510E	34–42	45	45	10	7–32	6.5	24
TGM454515E	34–42	45	45	15	7–32	6.5	24
TGM454520E	34–42	45	45	20	7–32	6.5	24

GORE® TAG® Conformable Thoracic Stent Graft with ACTIVE CONTROL System


Intended range of aortic inner diameters (mm)

16–42 mm
aortic diameter range
as few as **5 sizes**



31 x 26 mm

intended to treat proximal aortic diameter range of 24–29 mm and distal aortic diameter range of 19.5–24 mm



26 x 21 mm

intended to treat proximal aortic diameter range of 19.5–24 mm and distal aortic diameter range of 16–19.5 mm

Aortic accessories

GORE® DrySeal Flex Introducer Sheath

Catalogue number	Sheath size (Fr)	Minimum sheath ID (mm)	Nominal sheath OD (mm)	Working length (cm)
DSF1033	10	3.3	4.0	33
DSF1045	10	3.3	4.0	45
DSF1065	10	3.3	4.0	65
DSF1233	12	4.0	4.7	33
DSF1245	12	4.0	4.7	45
DSF1265	12	4.0	4.7	65
DSF1433	14	4.7	5.3	33
DSF1465	14	4.7	5.3	65
DSF1533	15	5.0	5.6	33
DSF1633	16	5.3	6.1	33
DSF1665	16	5.3	6.1	65
DSF1833	18	6.0	6.7	33
DSF1865	18	6.0	6.7	65
DSF2033	20	6.7	7.5	33
DSF2065	20	6.7	7.5	65
DSF2233	22	7.3	8.2	33
DSF2265	22	7.3	8.2	65
DSF2433	24	8.0	8.8	33
DSF2465	24	8.0	8.8	65
DSF2633	26	8.7	9.5	33
DSF2665	26	8.7	9.5	65

Aortic accessories

GORE® Molding & Occlusion Balloon Catheter

Catalogue number	Size (Fr)	Balloon diameter (mm)	Catheter length (cm)
MOB37	10	10–37	90

GORE® Tri-Lobe Balloon Catheter

BCM1634 – 10 mm Balloon Catheter

Intended aortic diameter (mm)	Recommended inflation volume (ml)
16	4
20	6
24	8
28	10
32	13

BCL2645 – 16 mm Balloon Catheter

Intended aortic diameter (mm)	Recommended inflation volume (ml)
26	12
30	14
34	18
38	21
42	25

 Consult Instructions
for Use
eifu.goremedical.com

Refer to *Instructions for Use* at eifu.goremedical.com for a complete description of all applicable indications, warnings, precautions and contraindications for the markets where this product is available. Rx Only

Products listed may not be available in all markets.

GORE, *Together, improving life*, ACTIVE CONTROL, C3, EXCLUDER, TAG and designs are trademarks of W. L. Gore & Associates.

© 2021 W. L. Gore & Associates GmbH 21327693-EN OCTOBER 2021

W. L. Gore & Associates, Inc.
goremedical.com

Asia Pacific +65 67332882

Australia/New Zealand 1800 680 424

Europe 00800 6334 4673

United States Flagstaff, AZ 86004 800 437 8181 928 779 2771

